Fassisi® BoTick

#01-09005-16
#01-09010-16
#01-09025-16

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EN Instruction Manual

The rapid test for direct detection of Borrelia bacteria in ticks

The Borrelia antigen test (Fassisi BoTick) is used for the direct detection of Borrelia bacteria (B. garinii, B. afzelii and B. burgdorferi sensu stricto) in ticks, the pathogens that cause Lyme disease. Lyme disease is a serious illness transmitted by ticks. The pathogen responsible for Lyme disease, Borrelia burgdorferi sensu lato, is a gram-negative, spiral-shaped and motile bacterium (spirochete). Worldwide, at least 13 genospecies belong to the Borrelia burgdorferi sensu lato complex. The three most important species, B. burgdorferi sensu stricto (Bss), B. garinii (Bg) and B. afzelii (Ba) are also pathogenic to humans.

Borrelia bacteria are transmitted by ticks of the genus Ixodes ricinus (common wood tick). In Europe, the prevalence of Borrelia bacteria in ticks varies from region to region (5–35%). In regions with a high prevalence of ticks, animals can become infected several times a year.

The clinical picture of Lyme borreliosis in dogs often develops gradually and is characterised by non-specific symptoms such as apathy, swollen lymph nodes or intermittent fever. In horses, symptoms such as lameness, joint swelling, neurological and eye changes have been observed. Even with effective antibiotic therapy, complete elimination of the pathogen is difficult, especially at an advanced stage of infection. Early antibiotic treatment is therefore essential for successful therapy. Early detection of Borrelia antigens allows non-specific symptoms, such as joint problems, that are suggestive of borreliosis to be better classified in the overall diagnosis, enabling rapid and targeted therapy to be initiated.

Sample information For tick use only

Carefully remove the tick, ensuring that the head is also removed. Tick tweezers are recommended for this purpose. Between removal and testing, the tick can be stored in the enclosed tick tube.

Note for optimal test results:

Optimal test results depend on the intensive crushing of the tick, especially its intestine. The Borrelia bacteria to be detected are located in the tick's intestine. Testing fresh, engorged ticks in combination with optimal crushing yields optimal test results.

Test procedure

Open the aluminum pouch, remove the test cassette, place the test cassette on a flat surface and unscrew the bottle of reagent and place it aside.

Figure 1

Place the tick in the tick tube and add 3 drops of reagent.

Figure 2

Crush the tick using the plastic stick; the intestine must be broken up and the fluid must change colour.

Use both the pointed and blunt ends of the plastic stick for this.

Figure 3

Remove the plastic stick and draw up the tick reagent mixture with the pipette.

Note: If the tick is heavily engorged and very large, the crushing process may result in the tick reagent mixture becoming very thick. In this case, dilute the mixture with additional drops of reagent until the liquid can be easily picked up with the pipette.

Figure 4

Place three drops of the tick reagent mixture onto the sample area of the test cassette. Make sure that you do not include any tick parts. These could clog the sample area. The liquid will begin to rise up the test strip.

If the liquid does not rise after 30 seconds, take the pipette and press the tip onto the sample area to activate the flow, or add another drop of reagent to the sample area.

Practical note

- Crush the ticks carefully, as the Borrelia bacteria to be detected are located in the tick's intestine.
- Do not place any tick parts on the sample field.
- Only use the plastic stick for the purpose described here.
- Remove the plastic stick from the reagent immediately after crushing.

Test result

The result of the test can be read after 10 minutes.

Positive test result (+)

If the result is positive, two red lines will appear in the evaluation field of the test cassette. A red line in the T area (T) indicates a positive test result. Borrelia bacteria have been detected in the tick.

The second red line in the C area (C) is the control line, which indicates that the test has been performed correctly. The C line is not a reference line and may show a different line intensity than the T line.

The fewer Borrelia bacteria present in the tick, the weaker the T line. This means that even a very weak test line represents a positive test result.

Negative test result (-)


If the result is negative, only the control line (C) will be visible in the evaluation field of the test cassette; no Borrelia bacteria could be detected in the tick.

Invalid test result

If no control line is visible, the test is invalid and must be repeated.

Components of the test kit

Number of components of the test kit and expiration dates: compare information on the outer packaging of the present test kit.

Storage of the test kit

The test kit must be stored between 2–30°C.

Please note before use

Use a new test cassette for every individual test. Only for one-time usage. For veterinary use only. Use only the original test components provided in the Fassisi kit. Use the test cassette within 60 minutes after opening the aluminum pouch. The test cassette must be in a horizontal position on a smooth surface under while the test is performed. Note the amount of sample material needed. An incorrect number of drops or too small drops may lead to false results. Consider the test results as invalid after the read out time. Do not use the test after the expiration date on the aluminum pouch. Dispose of all contaminated materials properly. Disinfect the work area after the test execution.

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