Testperformance	
Comparison Test 2016 (against ELISA):	
Sensitivity	Specificity
Fassisi TetaCheck	95,83 %
Tetanus-Antikörpernachweis	92,30 %

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EN Instruction Manual

To detect the tetanus antibody concentration

A tetanus infection, which occurs worldwide, is in most cases a wound infection caused by a previous injury. Contaminants, usually in combination with foreign bodies, get under the skin. There, the wound infection is caused by the toxins produced by the bacterium Clostridium tetani. These enter the central nervous system via the peripheral nerve tracts and can ultimately lead to death through respiratory paralysis. The course of the disease is very agonising for the animal and is almost always fatal.

For this reason, tetanus immunisation should be an absolute must for all horses. After successful basic immunisation, the vaccination must be regularly boosted in accordance with the recommendation of the StIKo Vet in order to maintain vaccination protection.

As scientific findings have shown that vaccination protection can last well beyond the vaccination interval specified by the vaccine manufacturer if basic immunisation has been carried out correctly, it is possible to check the existing tetanus antibody level using a rapid test.

After consistent, complete basic immunisation, the decision to vaccinate can be made dependent on the result of the rapid test. With the Fassisi TetaCheck, a quick check is available directly on site to regularly monitor this protection.

The tetanus antibody tests can be used, among other things - as a regular vaccination check - in the hospital sector - for serious injuries directly on site - in the breeding sector - at sporting events

Sample information

Serum and plasma
Recommended sample material is a freshly collected serum or plasma to achieve the highest detection sensitivity. Obtain serum or plasma as soon as possible after blood collection. Clear and non-hemolyzed sample material prevents slight background staining.




Whole blood
A whole blood sample should be used as quickly as possible. Heparin blood and EDTA blood may also be used. Hemolyzed samples should not be used for testing.

Good to know: If a sample from whole blood is used, this has a lower sensitivity than serum or plasma samples. If the test result is negative, a repeat test should be carried out with a serum or plasma sample, as the antibody titre limit is set with this sample material.

Test procedure

Open the aluminum pouch, remove the test cassette, place the test cassette on a flat surface and unscrew the bottle of reagent and place it aside.

Good to know: In the reaction field, you will see a green/blue line in the test and control region before the test is used. This is used for quality control and is washed away by the sample fluid during the test.

Practical note

If the fluid does not run up the test strip after 60 seconds, add an additional drop of the reagent into the sample well.

Test result

The result of the test can be read after 10 minutes.

Positive test result (+)
If the test result is positive, two red lines are visible in the evaluation field of the test cassette. The upper line (control line) confirms that the test has run correctly; the lower line (test line) indicates a positive test result. The tetanus antibody concentration is sufficient and is >0.1 IU/ml.

Even a weak test line is to be considered a positive test result, but is within the limit of the required antibody concentration.

Negative test result (-)

A red line (control line) is visible in the upper area of the evaluation field; no test line can be recognised. The tetanus antibody concentration is too low.


Note: The C-line is not a reference line and may have a different line intensity than the T-Line.

Invalid test result

If no control line appears after the test is conducted, the test is invalid.

Components of the test kit

Number of components of the test kit and expiration dates: compare information on the outer packaging of the present test kit.

Storage of the test kit

The test kit must be stored between 2-30°C.

Please note before use

Use a new test cassette for every individual test. Only for one-time usage. For veterinary use only. Use only the original test components provided in the Fassisi kit. Use the test cassette within 60 minutes after opening the aluminum pouch. The test cassette must be in a horizontal position on a smooth surface under while the test is performed. Note the amount of sample material needed. An incorrect number of drops or too small drops may lead to false results. Consider the test results as invalid after the read out time. Do not use the test after the expiration date on the aluminum pouch. Dispose of all contaminated materials properly. Disinfect the work area after the test execution.

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