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The fast measurement of foal IgG in blood samples

Foal losses in Germany, for example, are reported to be between 2% and 12%, depending on the phase of life, management and breeding area. In some countries, these figures are significantly higher. Identifying at-risk foals and assessing postnatal disease symptoms is often not easy, since waiting until clinically apparent changes occur costs valuable time.

The healthy development of newborn foals depends crucially on an adequate supply of immunoglobulins (IgG) from the mother's colostrum. An inadequate transfer of colostrum IgG is one of the most important predisposing factors for infectious foal diseases. Studies have shown that foal infection risk increases when the IgG concentration falls below 800 mg/dl.

In the context of a veterinary examination of every newborn foal, the importance of a routine determination of the IgG level as a way to identify animals that are more susceptible to disease cannot be overstated. It ensures the early diagnosis of an immunodeficiency and the initiation of further measures.

Sample information

Serum and plasma

Recommended sample material is a freshly collected serum or plasma to achieve the highest detection sensitivity. Obtain serum or plasma as soon as possible after blood collection. Clear and non-hemolyzed sample material prevents slight background staining.

Whole blood

A whole blood sample should be used as quickly as possible. Heparin blood and EDTA blood may also be used. Hemolyzed samples should not be used for testing.

Note: Whole blood samples have a lower detection sensitivity. If the result with whole blood is unclear, the test should be repeated with a serum or plasma sample.

Test procedure

Open the aluminum pouch, remove the test cassette, place the test cassette on a flat surface and unscrew the bottle of reagent and place it aside.

Good to know: In the reaction field, you will see a green/blue line in the control region before the test is used. This is used for quality control and is washed away by the sample fluid during the test.

Figure 1

Using the marked 20 µl pipette, aspirate the sample material up to the 20 µl mark.


Figure 2

Now add this material to the sample well and allow the 20 µl to be drawn in. Avoid the formation of air bubbles.

Figure 3

Take the opened reagent vial and add 2 drops of reagent to the sample well. After a few seconds, the liquid begins to run over the test strip. Avoid the formation of air bubbles.

Practical note

If the fluid does not run up the test strip after 60 seconds, add an additional drop of the reagent into the sample well.

Test result

The result of the test can be read after 10 minutes.

(+) No test line:

≥ 800 mg/dl

No test line is visible. The IgG concentration is high. It can be assumed that passive immunity has been successfully transferred (> 8 mg/ml).

(+/-) Weak test line:

400-800 mg/dl

If a line appears in the test area that is only weakly pronounced, the IgG concentration can be assessed as very low (between 4 and 8 mg/ml).

(-) Visible test line:

< 400 mg/dl

If a test line appears in the test area, it can be assumed with a high degree of probability that the transmission of passive immunity has failed. The IgG concentration is too low (< 4 mg/ml).

Note: The C-line is not a reference line and may have a different line intensity than the T-line.

Invalid test result


If no control line appears after the test is conducted, the test is invalid.

Good to know:

The more visible the test line, the lower the IgG concentration. If the test line is clearly visible, there has been an inadequate transfer of passive immunity. The IgG concentration is too low. Action is required.

Components of the test kit

Number of components of the test kit and expiration dates: compare information on the outer packaging of the present test kit.

Storage of the test kit

The test kit must be stored between 2-30°C.

Please note before use

Use a new test cassette for every individual test. Only for one-time usage. For veterinary use only. Use only the original test components provided in the Fassisi kit. Use the test cassette within 60 minutes after opening the aluminum pouch. The test cassette must be in a horizontal position on a smooth surface under while the test is performed. Note the amount of sample material needed. An incorrect number of drops or too small drops may lead to false results. Consider the test results as invalid after the read out time. Do not use the test after the expiration date on the aluminum pouch. Dispose of all contaminated materials properly. Disinfect the work area after the test execution.

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