Fassisi® BoDia Plus

#01-50305-16
 #01-50310-16
 #01-50325-16

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EN Instruction Manual

The advanced rapid test for calf diarrhoea

Diarrheal diseases have been the main cause of calf losses in the first weeks of life for many years. The spread of diarrhoea pathogens within a farm can be up to 100%. Losses of 30% are also not uncommon. The calves often suffer permanent damage, which manifests itself in the form of poorer performance or increased susceptibility to disease. In the case of diarrhoea, a distinction must be made between infectious and non-infectious causes. Triggers for infectious calf diarrhoea can be identified with antigen detection. If rota and coronaviruses are present, the intestinal mucosa is severely damaged and, in combination with bacterial pathogens such as E. coli K99, leads to increased water excretion via the intestine. Ultimately, this leads to a real "dehydration" of the calf. In addition to viruses and bacteria, cryptosporidia in particular are often a cause of severe diarrhoea.

However, not all pathogens have to be ingested by the calf via the environment. The spore-forming bacterium Clostridium perfringens occurs naturally in the intestinal flora. By creating imbalances in the intestinal flora, this bacterium can produce toxins, which can also cause diarrhoea in calves. These individual diarrhoea pathogens can cause severe diarrhoea alone or in combination.

The Fassisi BoDia Plus provides an overview of the current pathogen situation in the herd in just a few minutes; effective treatment can be initiated.

Sample information

The best sample material is fresh calf feces. The fecal sample must be at room temperature when the test is carried out.

Important note on sampling:
 – Do not take a fecal sample from the ground – Do not collect large particles
Hygiene note: Fecal samples can be infectious. Dispose of all contaminated materials in accordance with regulations and disinfect the work area after carrying out the test.

Sample procedure

Calf feces are sampled using a pipette, cotton bud or spoon, depending on the consistency of the feces.

For liquid feces: Use a pipette

Unscrew the cap of the sample tube. Pick up the liquid feces with the pipette. Add 3 to 4 drops of fecal material from the pipette into the sample tube.

The number of drops required depends on the consistency of the feces; the more liquid the feces, the more drops are required. Screw the sample tube closed again.

For creamy feces: Use a cotton swab

Unscrew the cap of the sample tube. Pick up the feces with the tip of the cotton swab. Hold the cotton swab in the sample tube and dissolve the sample material in the reagent by stirring. Remove the cotton swab and screw the sample tube closed again.

A lot does not help! Only the tip of the cotton swab should be covered with feces. If too much fecal material is used, the test strip will turn brownish. This can interfere with the sample run and lead to unclear or invalid test results.

For more solid feces: Use the spoon in the sample tube

Unscrew the entire upper part of the sample tube and remove it, including spoon. Pick up the feces with the spoon and smooth it out so that not too much fecal material gets into the sample tube.

Put the upper part of the sample tube back into the sample tube and screw it shut so that the sample material can dissolve in the reagent.

Test procedure

Before starting the testing procedure, open the aluminium packaging of the test cassette and place the cassette on a horizontal surface.

Note: In the evaluation field, there is a green line in the test and control region before the test is used. These are used for quality control and are washed away by the sample fluid when the test is used.

Figure 1

Take the sample tube with the fecal sample and shake it briefly.

Figure 2

Unscrew the cap of the sample tube.

Figure 3

Add three (3) drops of the sample material to each of the sample well of the test cassette. If the liquid does not run over the test strip after a few seconds, add another drop.

Practical note

In the event of severe running difficulties due to too much or slimy feces, the pipette can be pressed into the respective sample field to reactivate the test run.

Test result

The results of the test must be read after 10 minutes.

Positive test result (+)

The test line (T region) and control line (C region) become visible on the corresponding test strip. Even a weak test line is to be interpreted as a positive test result. In the figure on the left, rotaviruses were detected.

Negative test result (-)

Only the control line (C region) is visible on the corresponding test strip, the test is negative. No antigens could be detected in the sample. In the picture on the left, no coronaviruses were detected.

Invalid test result

No control line (C region) becomes visible. The test is invalid and should be repeated.

Interpretation aid

Weak test lines:

If a test line is visible, the test is positive. Even if the test line is only faintly visible.


Weak control lines:

Using too much feces may result in weak control lines.

There are some fecal samples that reduce the intensity of the control lines. Nevertheless, if a control line appears, the test is considered valid.

Components of the test kit

Number of components of the test kit and expiration dates: compare information on the outer packaging of the present test kit.

Storage of the test kit

The test kit must be stored between 2-30°C.

Please note before use

Only for single-use.
 For professional use only.
 Use a new test cassette and a new sample tube for each test.
 Only use the components supplied to carry out the test.
 Use the test cassette within one hour after opening the aluminum pouch.
 The test cassette must be in a horizontal position on a smooth surface under while the test is performed.
 Note the amount of sample needed. Too much feces may interfere with the run of the test. Please note the specified evaluation times.
 Do not use test cassettes after the expiry date.
 The feces may be infectious.
 Dispose of all contaminated materials in accordance with regulations and disinfect the work area after carrying out the test.

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