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


EN Instruction Manual

The bovine PAG pregnancy test on site

The Fassisi BoviPreg is a pregnancy test that detects the pregnancy-specific protein PAG (Pregnancy Associated Glycoprotein) in blood samples from cattle.

PAG is highly specific and is only produced after successful implantation of the egg in the animal's uterus. After it has been produced, the PAG enters the maternal bloodstream. This makes it possible to detect the PAG concentration in blood samples in the event of pregnancy. In contrast to progesterone, PAG is therefore a highly specific parameter for determining pregnant animals. The PAG concentration increases from day to day as the pregnancy progresses. The higher the PAG concentration in the animal's blood, the more reliable the pregnancy diagnosis. According to the literature, a PAG concentration > 2.0 ng/ml in serum from the 30th day after insemination is 94% certain to indicate pregnancy.

PAG concentrations > 2.0 ng/ml can be reliably detected using the Fassisi BoviPreg. The Fassisi BoviPreg offers an early pregnancy diagnosis without laboratory-based equipment, which is possible directly on the farm. This shortens calving intervals and optimises overall pregnancy management.

Sample information

Serum and plasma

Recommended sample material is a freshly collected serum or plasma (EDTA) to achieve the highest detection sensitivity. Obtain serum or plasma (EDTA) as soon as possible after blood collection. Clear and non-hemolyzed sample material prevents slight background staining.

Whole blood

Whole blood sample (EDTA) should be used as quickly as possible. Hemolyzed samples should not be used for testing.

Note: Whole blood samples have a lower detection sensitivity. In case of a negative test result with whole blood, the test should be repeated with a serum or plasma sample.

Test procedure

Open the aluminum pouch, remove the test cassette, place the test cassette on a flat surface and unscrew the bottle of reagent and place it aside.

Good to know: In the reaction field, you will see a green/blue line in the control region before the test is used. This is used for quality control and is washed away by the sample fluid during the test.

Figure 1

Using the marked 20 µl pipette, aspirate the sample material up to the 20 µl mark.


Figure 2

Now add this material to the sample well and allow the 20 µl to be drawn in. Avoid the formation of air bubbles.

Figure 3

Take the opened reagent vial and add 2 drops of reagent to the sample well. After a few seconds, the liquid begins to run over the test strip. Avoid the formation of air bubbles.

Practical note

If the fluid does not run up the test strip after 60 seconds, add an additional drop of the reagent into the sample well.

Test result

The result of the test can be read after 10 minutes.

Positive test result (+)

If the test result is positive, two red lines are visible on the test strip in the evaluation field of the test cassette. The upper line (control line) confirms that the test ran correctly; the lower line (test line) indicates the presence of PAGs.

Even a weak test line is to be regarded as positive evidence.

Note: Abortions can occur at any stage of pregnancy, especially in the first 40 days. An early pregnancy check must be confirmed by follow-up checks.

Negative test result (-)

Only a red line in the upper area of the reaction field (control line) becomes visible, no test line becomes visible. There are no PAGs in the sample at the current time or the PAG concentration is still too low. A repeat test of the PAG concentration after a few days is recommended as a control, as the PAG concentration increases over time in the event of pregnancy.

Note: The C-line is not a reference line and may have a different line intensity than the T-Line.

Invalid test result


If no control line appears after the test is conducted, the test is invalid.

Note regarding PAG concentration

If the PAG concentration in the serum is below 1.0 ng/ml, pregnancy can be ruled out with a high degree of probability and the rapid test is negative. At 1.0 to 2.0 ng/ml, no clear statement can be made and it is recommended to test another sample a few days later. A value of over 2.0 ng/ml, on the other hand, indicates a pregnancy with a high degree of certainty; the rapid test is positive.

Components of the test kit

Number of components of the test kit and expiration dates: compare information on the outer packaging of the present test kit.

Storage of the test kit

The Fassisi test kit may be kept and stored between 2-30°C. No cooling is required.

Please note before use

Use a new test cassette for every individual test. Only for one-time use. For veterinary use only. Use only the original test components provided in the Fassisi kit. Use the test cassette within 60 minutes after opening the aluminum pouch. The test cassette must be in a horizontal position on a smooth surface under while the test is performed. Note the amount of sample material needed. An incorrect number of drops or too small drops may lead to false results. Consider the test results as invalid after the read out time. Do not use the test after the expiration date on the aluminum pouch. Dispose of all contaminated materials properly. Disinfect the work area after the test execution.

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