

SKIN MARKER WITH RULER



Product

Marking the patient's skin is important in a range of clinical settings and for different reasons but the principal functions can be to:

- ensure surgery is performed on the correct site
- target non- invasive treatments
- plan the surgical approach
- ensure natural realignment of skin flaps during wound closure

MARKER



- ✓ Double head marker with regular and fine tip
- ✓ Large ink reservoir for longer product life
- ✓ Constant ink flow for clear marking
- ✓ Clip cap for regular tip and standard cap for fine tip

Gentian violet ink

- ✓ Colour works on all ranges of pigmented skin
- ✓ Quick drying, indelible, biocompatible
- ✓ Non-irritating, non-toxic
- ✓ Non-smearing, does not run



For quality
&
safety





RULER

Flexible 15 cm
In centimetres and inches
Ensures measurement even over contoured surfaces

Features

- Single-use, individually packaged in sterile peel pouches
- Latex-free
- Tapered tip for ergonomic handling

SMI tip

For best results, remove oil from the skin using alcohol before marking

Technical description

Package of 36
14.1 cm long (with caps)
Weight 7 g
High Density Polyethylene (HDPE)

Product	Unit	Order no
Skin marker sterile, with ruler	36	ZSM20

